(**DE**)

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(**EN**)

**Intended use**

A-silicone for the production of earmolds

**Indication for use**

For manufacturing of individual earmolds for hearing aids and hearing protection.

**Contraindication**

The product may not be used if the patient is allergic to one of the ingredients. There are no side effects to be expected from this medical device if processed and used properly. However, immune reactions (e. g. allergies) or local disorders of sensation (e. g. irritation in the auditory canal) cannot be ruled out in principle. If you become aware of adverse side effects, also in cases of doubt, please be sure to advise us in any case and describe the circumstances and symptoms as accurately as possible. We will investigate any notice.

**Target Patient Group**

Patients who need ear care.

**Intended users**

Technical professionals/earmold technicians

**Product description**

**Biopor® AB** is an addition curing two-component silicone with a final hardness of 25, 40 or 60 Shore A. **Biopor® AB light** with a final hardness of 25 Shore A and **Biopor® AB fluoreszent** with a final hardness of 40 Shore A are floatable, addition curing two-component silicones. These materials are suitable for the production of individual earmolds.

**Biopor® AB Xtrack** is an addition curing detectable 2-component silicone with a final hardness of 40 Shore A. It can be used for the production of detectable hearing protection earmolds, in particular for the food industry. Biopor® AB Xtrack is detectable by metal detectors and has the same effect as an iron ball with a diameter of 5 mm. It is also detectable by X-ray.

**Warning:** Metal detectors (or detection systems) have different default settings, depending on their learning process. Therefore, make sure that the Biopor® AB Xtrack material is detected by metal detectors (or detection systems) before first use.

**Processing**

Biopor® AB is automatically mixed from the double cartridge in a ratio of 1:1. For this purpose, insert the double cartridge into the injector. Remove the cap and dispense a little material from the double cartridge to ensure that the two components come out evenly. Then place a mixing canula onto the double cartridge and start filling the negative mold. Biopor® AB materials cure at room temperature. For easier removal from the FotoCast® form, we recommend using the separator before filling. To produce bubble-free earmolds, curing should be performed under pressure (4 bar) in the Polymax.

**Attention:** Biopor® AB light may only be used together with a separator in the FotoCast® process.

Single steps of the process are explained in "Production methods".

Please find processing times in "technical data" overview.

**Post-processing**

Take the earmolds out of the negative form and process them with cutting tools for soft materials. The earmolds can then be lacquered. To upgrade the surfaces, we recommend silicone lacquers, e. g. Lacquer B.eco.

**Disinfecting**

For thorough cleaning and disinfection, we recommend our OtoVita® care series.

**Cautions**

- Biopor® AB materials must not be used as ear impression silicone.
- The product is exclusively intended for the above-mentioned range of use and shall only be used by professionally qualified and instructed persons.
- Due to technical filling reasons there might be an unequal transportation of the two components of new cartridges. For this reason it is recommended to dispose the initially mixed pea-sized amount.
- Pre-convey the material before attaching the mixing canula.
- Polymerization, referring to the hardening of addition curing silicones, may be inhibited by latex gloves, creams, plastics, cleaning agents, etc. For this reason, we recommend that you wear commercially available polyethylene gloves.
- The setting times can be increased by lower temperatures or reduced by higher temperatures.
- Wear protective clothing! Set material is chemically resistant. Marks on your clothing cannot be removed.
- To avoid a cartridge cracking the material should be applied very slowly and with care when using a canula tip.
- Do not use canula tips with Biopor® AB (60 Shore A).
- Wear suitable personal protective equipment (protective gloves, safety goggles, face mask) during reworking of the cured material.
- After use, leave the mixing canula on the double cartridge as a cap.
- Never dispense into a mixing canula which has hardened.

**Keep away from children!**

**Note for one-time use**

Due to the physical features, the cured earmold silicone cannot be used several times.

**Disposal**

Disposal of contents/container in accordance with local/regional/national/international regulations and safety data sheet.

**Hazard identification**

Hazard and safety notes can be found on the product label and the corresponding safety data sheet.

**Serious incidents**

All serious incidents occurring in relation to the product shall be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is settled.

**Batch number / use-by date**

The batch number and the use-by date are located on both the outer packaging and the product. When submitting a complaint about a product, please always quote the batch number. Do not use the product after expiration of the use-by date.

**Dreve Otoplastik GmbH is not liable for any damages caused by improper application of the product.**

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Dreve Otoplastik GmbH Max-Planck-Straße 31 59423 Unna/Germany www.dreve.com