# Dubliersilikon transparent Dubliersilikon rosa

# Duplicating silicone transparent Duplicating silicone pink

egger Otoplastik +  
Labortechnik GmbH  
Aybühlweg 59  
87439 Kempten / Germany

egger

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## General product information

Addition-vulcanizing casting silicone for manufacturing negative molds in the production of earmolds and for the control molds.

**Duplicating silicone, transparent** is used primarily for processing UV-polymerizates or for producing control molds.

Final degree of hardness approx. 19–21 shore.

**Duplicating silicone, pink** is used in the acrylic technique. The components A and B can be mixed manually or mechanically. The negative mold is shrinkproof.

Final degree of hardness approx. 12 shore.


### Delivery units:

<b>Duplicating silicone, transparent</b>		<b>Duplicating silicone, pink</b>	
2 x 1 liter	REF 41000	2 x 1 liter	REF 41801
2 x 5 liter	REF 41001	2 x 5 liter	REF 41802
8 x double cartridges	REF 41003		

**Stability / Batch:** See tub / cartridge

### Storage:

Protect against: frost, heat, UV-radiation / sunlight. Keep container tightly closed in a cool, well-ventilated place.

### Indications:

<b>Duplicating silicone, transparent</b>	<b>Duplicating silicone, pink</b>
Component A – catalyst	Component A white – catalyst
Component B – base	Component B pink – base

- Avoid any contact between the components A and B before mixing them.
- Close the receptacle tightly after manual use.
- Avoid direct skin contact with the liquid material! Before handling the material protect your hands by wearing suitable gloves.
- Cross-linked materials are chemically resistant – avoid stains on your clothes.

### Notes for processing:

Before processing the material it must have reached a temperature of  $23\text{ °C} \pm 2\text{ °C}$  ( $73.4\text{ °C} \pm 35.6\text{ °F}$ ) at a relative air humidity of  $50\% \pm 5\%$ .

Processing time:	4–6 minutes
Curing time without pressure pot:	15–30 minutes
Curing time with pressure pot:	10–15 minutes at $40\text{ °C}$ ( $104\text{ °F}$ ) and 2–5 bar

Lower temperatures prolong and higher temperatures reduce time of processing and curing.

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






## Processing

### Manual mixing

Mix the two components A and B (mixing ratio 1:1) until material has homogeneous color. Approx. 40–50 g of duplicating silicone is needed for one negative mold.

### Automatic mixing

The material can be processed with the egger silicone-mixing device ESG 2 (REF 54560) or with an inject pistol (REF 25100). Mix the components A and B bubble-free in the mixing canula (mixing ratio 1:1). Do not use mixing canula (REF 25153 or REF 61402ff) if it contains hardened material.

### Filling

When filling the material into the investment mold make sure that the tip of the mixing canula always stays in the material to avoid the formation of bubbles.

### Polymerization

Afterwards let the material cure at room temperature or put the filled investment mold in the pressure pot (REF 54200/49800) for 10–15 minutes at 2–5 bar for a bubble-free polymerization. After the polymerization of the duplicating silicone, use the impression removing instrument (REF 46700) to carefully demold the impression. Before continuing with processing the impression, check the negative mold for wax remainders and remove them.

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7908607-2014-12

**Please follow the product information  
and the material safety data sheet.**