Service- und Reparaturset egger eLED.LP3

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Wendemotor (REF 89320)
Serien-Nr.:

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oder Wendemotor (REF 89330)
Serien-Nr.:

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Service and Repair Kit egger eLED.LP3

Operating instructions

REF 31700 (100 -240 V)

REF 31710 (100 -240 V)

Service and Repair Kit egger eLED.LP3

Please read the operating instructions carefully before working with the device!

Please keep the operating instructions in a safe place for future reference!

Index	page
1 General	11
2 Scope of delivery, Accessories, Technical Data	11
2.1 Scope of delivery of the eLED.LP3	12
2.2 Accessories	12
2.3 Technical data of the eLED.LP3	12
3 Safety Regulations	12
3.1 For a Safety Conscious Work Environment	12
4 Installation, Handling, Initial Operation	13
4.1 Assembling and Connecting the Unit	13
4.2 Initial Operating and Handling	13
4.3 Advice for Disposal	13
5 Working with UV-Polymerisates	14
5.1 Specific Storage Regulations	14
5.2 Specific Processing Regulations	14
5.3 Lacquering with LP/H laquer, LP/H laquer antibac	14
5.4 Lacquering with LP/W lacquer	14
5.5 Modifications/Repairs with Modeling Pastes	15
5.6 Optional manufacturing of BTE earmolds and ITE shells	15
5.7 Processing Times of LP/H Materials	16
6 Warranty/Service/After-Sales Service	16

Caption

Important Information!

Warning hot surfaces – do not touch the LED reflectors – risk of burns!

Warning UV radiation – do not look into the radiation source – risk of damage to the eyes!

Class II equipment.

Caution!

Service and Repair Kit egger eLED.LP3

1 General

The egger service- and repair kit is suited for a number of uses in the acoustician shop, during home visits and in the otoplastics laboratory. Service operations can be executed immediately right in the shop – fast, efficiently, no lengthy preparation needed. So you can enjoy more independence while offering your customer the best possible service.

All possible uses and advantages at one glance:

Modeling pastes

- Repairs of BTE and ITE earmolds that are made of UV-polymerizates or other hard acrylates
- Minimizing or sealing ventilation bores
- Post-production attachment of retention points or a lip seal
- Modifications of vents on ITE earmolds
- Handy 20 g tubes for efficient and clean processing

Lacquers

- Coatings with LP/H lacquer or LP/H lacquer antibac
- Post-production varnishing (“glazing”) of polished earmolds using LP/H lacquers – a preventive measure against skin irritation
- Partial build-up of earmolds to improve their seal
- Adhesive coating (“non-slip coating”) using LP/W lacquer to improve an earmold’s hold
- Usable for all UV-polymerizates and hard acrylates
- Available in handy 20 ml or 50 ml bottles including lacquering brush

The light-curing unit egger eLED.LP3 has been especially developed for the above mentioned applications and is characterized by

- Maintenance-free working
- Durable LED light sources
- Easy handling
- Simple cleaning and care of the powder-coated metal casing

This egger service and repair kit eLED.LP3 is rounded off by the necessary accessories for a quick gluing in of sound tubes.

Furthermore, we can offer a variety of devices and useful accessories from our wide range of products such as:

- modeling pastes, lacquers and accessories
- cutters, drills and tools
- technical units
- compressors for compressed air
- compressed-air sprays

2 Scope of delivery, Accessories, Technical Data

2.1 Scope of Delivery

Please check whether the delivery is complete and whether all parts are free from transport damage. Immediately inform us or the distributor about missing parts or transport damage.

The service and repair kit eLED.LP3 consists of:

1 light unit egger eLED.LP3	REF 51350
1 turning motor battery-operated (kit 31700)	REF 89320
or 1 electrical turning motor (kit 31710)	REF 89330
1 tube LP/H modeling paste transparent 20 g	REF 30800
1 tube LP/H modeling paste beige 20 g	REF 31000
1 tube LP/H modeling paste red 20 g	REF 29610
1 tube LP/H modeling paste blue 20 g	REF 29615

Service and Repair Kit egger eLED.LP3

1 bottle LP/H lacquer antibac 50 ml	REF 30621
1 bottle LP/W lacquer 50 ml	REF 30700
1 bottle LP-cleaner 100 ml	REF 31301
1 bottle surface cleaner 250 ml	REF 31600
1 hand-spray pump for surface cleaner	REF 86732
1 jar round, for LP-cleaner	REF 52800
1 meter nylon tube as support for lacquering	REF 97200
1 modeling instrument	REF 31215
1 timer to control polymerization times	REF 68601
1 Uniglu instant adhesive	REF 39501
5 measuring tips large, for Uniglu	REF 39504
5 measuring tips small, for Uniglu	REF 39505
10 sound tubes pre-bent, 2 x 3.1 mm	REF 95100
1 user and processing instruction	

2.2 Accessories

Insert and metalized carton	REF 50351
Metalized carton (5 units) with adhesive dots	REF 51352
Protective film, transparent, self-adhesive	REF 51353

2.2 Technical Data o

Power supply: $U = 100 - 240 \text{ V (AC)}$

Mains frequency: $f = 50/60 \text{ Hz}$

Power input: $P = 22 \text{ VA}$

Life expectancy of the LED bulbs: 20.000 operating hours

Connection: Euro plug

3 Safety Regulations

The electrical light-curing unit egger eLED.LP3 is used for the polymerization of UV-curing modeling pastes and lacquers. A use other than described in the operating instructions may lead to injuries of the operator or of third persons.

Furthermore, the device or other material assets may be damaged. The device must be used in technically perfect condition only and the safety regulations must be strictly complied with.

Warning hot surfaces – do not touch the LED reflectors – risk of burns!

Warning UV radiation – do not look into the radiation source – risk of damage to the eyes!

3.1 For a Safety Conscious Work Environment

- Protect the unit from dust. Use the unit on dry working places only.
- Store the unit in closed rooms.
- Qualified and trained persons must use the device only.
- Children must be kept away from the unit.
- Surrounding temperature between $10 \text{ }^\circ\text{C}$ and $30 \text{ }^\circ\text{C}$. Relative air humidity between 10 % and 70 %.
- Before putting the device into operation, the user must read the operating instructions and ask the manufacturer in case of any questions.
- Use the unit in perfect and clean condition only and by observing the operating instructions.
- Connect the device to permissible sockets only.
- Use the included cable only.
- Please adhere to the acknowledged safety regulations.
- This device is used for the processing of UV-curing materials in a spectral range of 360 – 400 nm.

Service and Repair Kit egger eLED.LP3

- Please also adhere to the processing instructions for each material used.
- Disconnect the unit from electricity supply prior to any service or maintenance.
- Do not pull the mains plug out of the socket at the mains cable.
- Should the mains cable be damaged, it must be replaced at once. Defective mains cables must not be used.
- The unit may not be opened by operators. There are no interchangeable parts such as fuses inside.
- Repairs may only be executed by the manufacturer.

Advice for cleaning: Clean the unit with a moist cloth only – do not use solvents, do not douse the unit into water.

On self-accomplished repairs or if other spare parts are used, all warranty and liability claims will become extinct!

4 Installation, Handling, Initial Operation

4.1 Assembling and Connecting the Unit

The electrical light-curing unit egger eLED.LP3 corresponds to the state-of-the-art standards and to the acknowledged safety regulations. Nevertheless dangers can arise. Please observe the following, before putting the unit into operation:

- The device is to be positioned in a way that a fast separation of mains is possible very easily and at any time.
- Position the device horizontally and place it in its suitable working position.
- Protect the unit against vibration.
- Check whether the unit is dry.
- Check whether the device is in technically perfect condition, especially the mains cable!
- Plug in the mains cable.

4.2 Initial Operating and Handling of the Light-Curing Unit egger eLED.LP3

Set the timer to the desired time (1, 3 or 5 minutes) by pressing the timer button and then start the polymerization process by pressing the start button.

Meaning of illuminating diodes

- 1 minute •
- 2 minutes • •
- 5 minutes • • •

Stopping the polymerization process

Keeping the start-stop-button pressed longer, the polymerization process can be stopped. In case of not using the device for 30 minutes, it changes into state of rest. The illuminating diodes turn off.

4.3 Advice for Disposal

The symbol on the unit means, that this unit may not be disposed in general trash, but must be disposed according to regional recycling guidelines for electronic devices and applicable environmental and disposal rules.

For further information on disposal and recycling please contact your regional authorities or your distributor of the unit.

Service and Repair Kit egger eLED.LP3

5 Working with the UV-Polymerizates

5.1 Specific Storage Regulations

The materials are to be stored in a cool and dry place (see processing instructions for the various materials). Please mind that the containers are closed at any time, when not in use and that they are closed immediately after use. Otherwise, there might be danger of activating an unintentional chemical reaction through the permeability of light (curing of the material). The expiry date, applicable when stored adequately, is indicated on the product.

TIP: We recommend keeping already opened packages of lacquers or instant adhesives in the refrigerator.

5.2 Specific Processing Indications

Please avoid direct skin contact with the materials. Before using them protect your hands by using special skin protection products or wearing gloves.

The materials are to be used by authorized and trained staff only. Processing of the lacquers should be done underneath a suction unit. Always use the lacquering brushes for one type of lacquer only.

Details on the disposal for each product can be taken from our material safety data sheets, which can be downloaded in PDF format from our website:

www.egger-labor.com.

5.3 Lacquering with LP/H Lacquer – LP/H Lacquer antibac

For the surface finish all at present known hard acrylate-based otoplastic materials can be coated with LP/H lacquer.

Using the LP/H lacquer antibac, you can protect hard earmolds with a special antibacterial coating with high gloss effect.

- The entire surface of the earmold always has to be roughened with corundum paper, grit 180, (REF 40704) first.
- Then, fix the BTE earmold onto a piece of nylon tube (REF 97200). For an ITE earmold, the alligator clip of the turning motor can be fixed on the opened battery door.
- Then, clean and degrease the surface carefully with surface cleaner by using the hand-spray pump (REF 86732).
- Apply the lacquer evenly with the brush. Alternatively, the earmold can also be immersed in a dipping jar (REF 31400 or 31401) filled with lacquer; but do not immerse ITE earmolds with built-in electronics or face-plates! Then, make sure that excess lacquer is shaken off.
- Let the lacquer dry in the light-curing unit rotating for about 2 – 3 minutes. Afterwards the earmold is directly ready for use.

5.4 Lacquering with LP/W Lacquer

For the surface finish all at present known hard acrylate-based otoplastic materials can be coated with LP/W lacquer to improve the fitting of the earmold (“non-skid coating”).

- The entire surface of the earmold always has to be roughened with corundum paper, grit 180, (REF 40704) first.

Service and Repair Kit egger eLED.LP3

- Then, fix the BTE/ITE earmold onto a piece of nylon tube (REF 97200). For an ITE earmold, the alligator clip of the turning motor can be fixed on the opened battery door.
 - Then, clean and degrease the surface carefully with surface cleaner by using the hand-spray pump (REF 86732). Apply the lacquer evenly with the brush.
 - Let the lacquer dry in the light-curing unit rotating for about 6 – 9 minutes.
 - After the polymerization, the earmold (but do not immerse ITE earmolds with built-in electronics) is put into the LP-cleaner bath (in the round jar REF 52800) and then a second polymerization in the light-curing unit egger eLED.LP3 is to be done for about 6 minutes. This process is necessary to remove the inhibition layer (sticky surface) as a result from the first polymerization.
 - To finish the process, the earmold is rinsed under water and is blown dry.
 - Then, clean earmolds for 3 minutes in ultrasonic bath of water / isopropanol 1:1 (weighting %), treat with isopropanol and let dry at 37 °C for approx. 60 minutes.
- Then, depending on the coating thickness applied, let it cure for about 2 – 6 minutes in the light-curing unit egger eLED.LP3.
 - After the polymerization of the surface, wipe off the earmold by using a surface cleaner moistened cloth.
 - If required, rework the earmold with a cutter and prepare it for lacquering as described in 5.3 / 5.4

5.6 Optional Manufacturing of BTE Earmolds and ITE Shells

Using the optionally available metalized insert, processing of egger LP/H UV-polymerizates for the manufacturing of BTE earmolds and ITE shells is possible.

Tip: the metalized insert carton is fixed with included adhesive dots at insert. For this purpose, adhesive dots are attached to punched marking of metalized carton and pressed to the insert.

5.5 Modifications / Repairs with Modeling Pastes

The LP/H modeling pastes can be universally used for various modifications and repairs with all kinds of hard UV-curing polymerizates or acrylates.

- First of all, the part of the earmold that is to be processed always has to be roughened with corundum paper, grit 180, (REF 40704).
- Then, clean and degrease the surface carefully with surface cleaner by using the hand-spray pump (REF 86732).
- Apply the modeling paste directly from the tube with the modeling instrument (REF 31215) and smooth the transitions to the ground.

Service and Repair Kit egger eLED.LP3

5.7 Processing Times of the Individual LP/H Materials

LP/H Material	1st Polymerization	2nd Polymerization	3rd Polymerization
LP/H transparent	2 min. + black foil	3 min.	5 min. LP Cleaner
LP/H pink-transparent	2 min. + black foil	3 min.	5 min. LP Cleaner
LP/H beige	1 min. + black foil	1 min.	5 min. LP Cleaner
LP/H flesh-colored	1 min. + black foil	1 min.	5 min. LP Cleaner
LP/H blue	15 sec. + black foil	1 min.	5 min. LP Cleaner
LP/H red	15 sec. + black foil	1 min.	5 min. LP Cleaner

6 Warranty / Service / After-Sales Service

The warranty period for our devices starts with the date of delivery and is 12 months. Excluded from the warranty are modifications of the unit such as lamp modules, reflectors as well as all other parts of the unit which are subject to wear and tear.

In case of improper handling of the electric light-curing unit all warranty and liability claims will become extinct!

We will be pleased to be of assistance to you in case of damage or other malfunctions:

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E-mail international: sales@egger-labor.de

Internet: www.egger-labor.com

Please indicate the following data:

Device type:

Electric light-curing unit egger eLED.LP3 (REF 51350)

serial no.:

.....
(see label)

Turning motor (REF 89320)

serial no.:

.....
(see label)

or Turning motor (REF 89330)

serial no.:

.....
(see label)

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egger Otoplastik + Labortechnik GmbH