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Intended use
Light-curing lacquer for earmolds

Extended purpose
Light-curing lacquer for hard and soft earmolds

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Indication for use
Thermosoft Lack is a light curing lacquer for coating BTE earmolds and ITE shells.

Contraindication
If an acrylic allergy exists, an allergic reaction may occur.

Patient target group
People with impaired hearing ability to wear a hearing aid.

Intended users
Qualified personnel from earmold laboratories and hearing aid acousticians

Requirements / equipment
Polylux 1000N (without nitrogen), Gysol, Oxystop Fluid

Basic Material Properties
See technical datasheet (TDS). Available on request.

Specific Manufacturing Considerations
Caution
During the application we recommend wearing personal protective equipment i.e. suitable gloves, safety goggles etc. Use a suitable extraction system or ensure sufficient ventilation when working with the lacquer.

Further processing
Preparation
The surfaces of the BTE earmold or the ITE shell should be free from dust and grease residues. Roughening of the plastic surface is recommended for improved wetting behaviour.

Further processing
Thermosoft Lack is a lacquer for surface treatment of soft earmolds (Thermosoft) and hard earmolds (FotoTec® or Fotoplast®) which remains soft after polymerisation. The lacquer can be applied using a dipping method or with a brush. Dip the BTE earmolds or ITE shells to be lacquered into the

Lot number / Expiry date
The batch number and the use-by date are located on both the outer packaging and the product. When submitting a complaint about a product, please always quote the batch number. Do not use the product after expiration of the use-by date.

Thermosoft Lack and allow to drip dry. Take care to ensure that no air bubbles form. The lacquer film should be as thin as possible, however covering.

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