# egger EL 4

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Operating Instruction

REF 51300

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English Version Page 10 – 16

# egger EL 4

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## Operating instructions

egger EL 4  
REF 51205 (230 V)

Please read the operating instructions carefully before working with the device!

Please keep the operating instructions in a safe place for future reference!

Index	page
<b>1 General</b>	<b>11</b>
<b>2 Safety regulations</b>	<b>11</b>
2.1 For a safety conscious work environment	12
<b>3 Installation, handling and initial operation</b>	<b>12</b>
3.1 Scope of delivery	12
3.2 Technical data	12
3.3 Spare parts and accessories	12
3.4 Assembling and connecting the unit	12
3.5 Starting-up and operation	13
<b>4 Working with UV-polymerizates</b>	<b>13</b>
4.1 Specific storage regulations	13
4.2 Specific processing notes	13
4.3 Lacquering with LP/H lacquer antibac	13
4.4 Lacquering with LP/W lacquer	14
4.5 Modifications/repairs with modeling pastes	14
<b>5 Maintenance</b>	<b>15</b>
5.1 Exchanging light tubes	15
5.2 Work performance of the light tubes	15
5.3 Exchanging of reflecting foil	15
<b>6 Warranty/Service/After-sales service</b>	<b>15</b>

## Caption

Important Information

Caution

## 1 General

The proven egger light unit EL 4 for the universal use in the otoplastic laboratory. It was especially developed to polymerize UV-curing resins, lacquers and modeling pastes for lacquering (vitrification), modification or for repairs.

The size of the polymerization chamber is designed so that up to 25 earmolds or ITE shells can be processed in only one work step.

The egger EL 4 is exclusively equipped with UVA light tubes to achieve an extra high emission maximum for the polymerization surface. Due to the number and order of the six UVA light tubes a high efficiency level is reached. A specific fan prevents the light chamber from over-heating.

NEW: The current polymerization times of two and four minutes are preset. For this reason the manual control of the polymerization times is omitted and the manufacturing processes are becoming even more comfortable. The transition button “2 min”/“4 min” is used to preselect the polymerization times. The light tubes automatically switch themselves off after termination of the selected polymerization time.

The light-curing unit egger EL 4 was especially developed for the applications described below and is characterized by

- Maintenance-free working
- Fast exchanging of the UV tubes
- Easy handling
- Easy to clean powder-coated metal casing

In the following a brief overview about the applications and the materials available:

### LP/H Modeling pastes

- Repairs of BTE and ITE earmolds made of UV-polymerizates or other hard acrylates
- Minimizing or sealing ventilation bores
- Post-production attachment of retention points or a lip seal
- Modifications of vents on ITE earmolds

### Lacquers

- Applying of an antibacterial coating with LP/H lacquer antibac
- Post-production varnishing (“glazing”) of polished earmolds using LP/H lacquers – a preventive measure against skin irritation
- Partial build-up of earmolds to improve their seal
- Adhesive coating (“non-slip coating”) using LP/W lacquers to improve an earmold’s hold
- Usable for all UV-polymerizates and hard acrylates

Furthermore, we can offer you a variety of devices and usable accessories of our wide range of products. Ask for our catalog or visit us on our homepage [www.egger-labor.com](http://www.egger-labor.com).

## 2 Safety regulations

 The electrical light-curing unit egger EL 4 is used for the polymerization of UV-curing modeling pastes and lacquers. A use other than described in the operating instructions may lead to injuries of the operator or of third persons. Furthermore, the device or other material assets may be damaged. Therefore, the device must be used in technically perfect condition only and the safety regulations must be strictly complied with.

# egger EL 4

## 2.1 For a safety conscious work environment

- Protect the unit from dust. Use the unit on dry working places only.
- Qualified and trained persons must use the device only.
- Before putting the device into operation, the user has to read the operating instructions and ask the manufacturer in case of any questions.
- Use the unit in perfect and clean conditions only and with full operational light tubes.
- Connect the device to permissible sockets only.
- Please adhere to the safety regulations accepted.
- This device is used for the processing of UV-curing materials in a spectral range of 360 - 390 nm. Please also adhere to the processing instructions and product information for each material being used.
- Do not pull the mains plug out off the socket at the mains cable.
- Should the mains cable be damaged, it must be replaced at once. Defective mains cables must not be used.
- The device must be repaired by the manufacturer only.
- The unit must be switched off during service and maintenance operations (cleaning, replacing of UV light tubes or reflecting foil). Pull the mains cable out of the mains plug!
- On replacing use the original egger light tubes and reflectors only.
- Do not use any dissolvent when cleaning the device. Do not immerse the unit in liquids.
- Children must be kept away from the unit!

 On self-accomplished repairs or if other spare parts are used, all guarantee and liability claims will become extinct!

## 3 Installation, handling and initial operation

### 3.1 Scope of delivery

Please check whether the delivery is complete and whether all parts are free from transport damage. Immediately inform us or the distributor about missing parts or transport damage.

The UV-curing unit egger EL 4 consists of:

1 varnishing drawer

1 mains cable

1 operating and processing instructions

and therefore can be put into operation right away.

### 3.2 Technical data of the light-curing unit egger EL 4

Power supply: 230 V +/- 10%

Mains frequency: AC 50 Hz

Connection: mains plug Schuko (German)

Power input: AC 120 Hz

### 3.3 Spare parts and accessories for the light-curing unit egger EL 4

UVA – light bulb REF 50103

Reflecting foil set REF 51209

Varnishing drawer REF 51208

### 3.4 Assembling and connecting the unit

The electrical light-curing unit egger EL 4 corresponds to the state-of-the-art standards and to the acknowledged safety regulations, nevertheless dangers can arise. Please observe the following, before putting the unit into operation:

- Qualified and trained persons must use the device only.
- Please adhere to the safety regulations accepted.
- Use the unit in perfect and clean conditions only and by observing the operating instructions.
- Position the device horizontally and place it in its suitable working position.

- Protect the unit against impact.
- Check whether the unit is dry.
- Check whether the device is in technically perfect condition, especially the mains cable!
- Plug in the mains cable.
- 

### 3.5 Starting-up and operation of the light-curing unit egger EL 4

All operating and control elements are fixed on top of the unit.

The mains switch (Power) is used to switch the unit on and off. On longer operating breaks the device should be switched off.

The transition button “2 min”/“4 min” is used to preselect the polymerization times. The light tubes automatically switch themselves off after termination of the selected polymerization time.

The START button switches on the UV light tubes and the fan. With the STOP button the light tubes can be switched off at any time.

The unit is designed for the operation with UVA light tubes. The light tubes are already mounted.

## 4 Working with the UV-curing materials

### 4.1 Specific storage regulations

The materials are to be stored in a cool and dry place (at a maximum room temperature of 23° C / 73.4° F). Please mind that the containers are closed at any time, when not in use and that they are closed immediately after use. Otherwise, there might be the danger of activating an unintentional chemical reaction through the permeability of light (curing of the material). The expiry date, applicable

when stored adequately, is indicated on the product.

**TIP:** We recommend keeping already opened packages of lacquers in the refrigerator.

### 4.2 Specific processing notes

Please avoid direct skin contact with the using them protect your hands by wearing gloves or using special skin protection products materials. Before.

The material is to be used by authorized and trained staff only. Processing of the lacquers should be done beneath an absorbing device (e.g. egger S/AB Timer REF 61400). Use the lacquering brushes always for one type of lacquer only. Please adhere to the processing notes for each material.

Details on the disposal for each product can be taken from our material safety data sheets, which can be downloaded in PDF format from our website: [www.egger-labor.com](http://www.egger-labor.com).

For processing, we recommend using additional accessories the egger lacquering holder kit REF 51207, consisting of magnetic lacquering holders (REF 51206), Erko-gum (REF 30911) and a nylon tube (REF 97200).

### 4.3 Lacquering with LP/H lacquer antibac

By using the LP/H lacquer antibac (REF 30621ff), you can protect hard earmolds with an anti-bacterial coating with high gloss effect

- First, the part of the otoplastic that is to be processed always has to be roughened with Corundum paper, grit 180, (REF 40704).

## egger EL 4

- Then, fix the BTE earmold onto a piece of nylon tube (REF 97200) and then onto the lacquering holder. For an ITE earmold, the special Ergokum can be fixed directly on the opened battery door.
- Then, clean and degrease the surface of the otoplastic carefully with Surface Cleaner (REF 31600ff).
- Apply the lacquer evenly. Alternatively, you can also immerse the otoplastic (but no ITE device with a built-in electronic or face-plates!) into a dipping jar filled with lacquer (REF 31400 or 31401). Make sure that spare lacquer is shaken off.
- Then let the lacquer cure in the light-curing unit for 2 minutes. Afterwards the otoplastic is directly ready to use.

### 4.4 Lacquering with LP/W lacquer

For the surface finish the light-curing soft otoplastics (made of LP/W) are always lacquered with LP/W lacquer (REF 30700ff). Furthermore, all at present known hard acrylate-based otoplastic materials can be coated with LP/W lacquer to improve the fitting of the earmold (“non-skid coating”).

- First, the part of the otoplastic that is to be processed always has to be roughened with Corundum paper, grit 180, (REF 40704).
  - Then, fix the BTE earmold onto a piece of nylon tube (REF 97200) and then onto the lacquering holder. For an ITE earmold, the special Ergokum can be fixed directly on the opened battery door.
  - Then, clean and degrease the surface of the otoplastic carefully with Surface Cleaner (REF 31600ff).
  - Apply the lacquer evenly and cure it for 2 minutes in the light-curing unit.
- After the polymerization, the otoplastic (but do not immerse ITE earmolds with built-in electronics – see note below!) is polymerized with LP-cleaner for about 4 minutes. This is necessary to remove the inhibition layer (sticky surface).
  - Then, the otoplastic is rinsed under water and blown dry.

ITE earmolds with a built-in electronic can also be coated with LP/W lacquer,. Please observe, however that the otoplastic should not be put into the LP-cleaner bath! To remove the inhibition layer use a Surface Cleaner moistened cloth and wipe off the surface thoroughly.

### 4.5 Modifications/repairs with modeling pastes

The LP/H modeling pastes can be universally used for various modifications and repairs with all kinds of hard UV-curing polymerizates or acrylates. These are available in the colors trans-parent, beige, red and blue.

- First, the part of the otoplastic that is to be processed always has to be roughened with Corundum paper, grit 180, (REF 40704).
- Then, clean and degrease the otoplastic with Surface Cleaner (REF 31600ff).
- Apply the modeling paste directly from the tube with the modeling instrument (REF 31215) and smooth the transitions to the ground.
- Then, depending on the coating thickness applied, let it cure for about 2 or 4 minutes in the light-curing unit.
- After the polymerization of the surface, wipe off the otoplastic by using a Surface Cleaner moistened cloth.
- If required, rework the otoplastic with a cutter and prepare it for lacquering as de-scribed in 4.3/4.4.

## 5 Maintenance

 Before each maintenance, the device must be unplugged from the mains supply!

### 5.1 Exchanging light tubes

On each exchanging of the tubes (REF 50103) always exchange them all. Please be sure always to have some spare tubes in stock.

- Pull the tubes evenly backwards out of the socket. Do not cant!
- Just put the new tubes into the socket until they have snapped into place.
- The unit is now ready for action and can be connected to the mains supply.
- Used light tubes are special waste and the legal provisions for their disposal have to be adhered to.

### 5.2 Work performance of the halogen lamp

The work performance of the light tubes depends decisively on the frequency of switching on the tubes and of the mains voltage. According to information from the manufacturer, the work performance averages approx. 1000 hours.

As the burning time of the light tube per polymerization is fairly short, we recommend:

- exchanging the tubes every ½ year, if the unit is in use for approx. 8 hours per day
- exchanging the tubes every year, if the unit is in use for not more than 4 hours per day

Abrasion of the light tubes can be noticed on elongated polymerization times and a dark discoloration along the plug and socket area. Used light tubes are hazardous waste. Please observe the current legal regulations.

### 5.3 Exchanging of reflecting foil

The polymerization performance of the light-curing unit depends on the cleanliness of the light chamber. The built-in reflectors (REF 51209) amplify the irradiation and always have to be clean, otherwise, the polymerization times can change. Therefore, exchange the foils in time.

## 6 Warranty/Service/After-sales service

### Warranty

The warranty period for our devices starts with the date of delivery and is 12 months. Excluded from the guarantee are modifications of the unit such as lamp modules, reflectors as well as all other parts of the unit which are subject to wear and tear.

In case of improper handling of the electric lightcuring unit egger EL 4 all warranty and liability claims will become extinct!

We will be pleased to be of assistance to you in case of damage or other malfunctions

egger  
 Otoplastik + Labortechnik GmbH  
 Aybühlweg 59  
 87439 Kempten/Germany  
 Phone international: +49 831 58113-60  
 Fax international: +49 831 58113-14  
 E-mail: sales@egger-labor.de  
 Internet: www.egger-labor.com

Please indicate the following data:

Device type:  
 Electric light-curing unit egger EL 4 (REF 51205)

.....  
 (see label)