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Light curing lacquer for earmolds

Extended purpose

Light curing lacquer for BTE earmolds and ITE shells

Indication

Hearing loss or impaired hearing

Indication for use

NanoScreen™ Lack is a light curing lacquer for coating BTE earmolds and ITE shells.

Contraindication

If an acrylic allergy exists, an allergic reaction may occur.

Patient target group

People with impaired hearing ability to wear a hearing aid.

Intended users

Qualified personnel from earmold laboratories and hearing aid acousticians

Requirements / equipment

Polylux 1000N (with nitrogen), Sonolux Flash Evo (with nitrogen), PCU LED N₂ (with nitrogen)

Basic Material Properties

See technical datasheet (TDS). Available on request.

Specific Manufacturing Considerations**Caution**

During the application we recommend wearing personal protective equipment i.e. suitable gloves, safety goggles etc. Use a suitable extraction system or ensure sufficient ventilation when working with the lacquer.

Process description**Preparation**

The surfaces of the BTE earmold or the ITE shell should be free from dust and grease residues. Roughening of the plastic surface is recommended for improved wetting behaviour.

Further processing

NanoScreen™ Lack is a light-curing single-component lacquer. After polymerisation, the result is a homogeneous and highly abrasion-resistant low-energy surface which reduces the adhesion of biofilms to the plastic surface of the earmold. Take care to ensure that no air bubbles form. For ITE shells

with integrated electronics, the lacquer should be applied in a very thin coating with a brush. The lacquer film should be as thin as possible. If you wish to achieve a thicker coating, briefly polymerize the first, thin lacquer film, then apply a second layer. Turn the earmold as required.

Post-curing

Polylux 1000N (with nitrogen)	2 min
Sonolux Flash Evo (with nitrogen)	1 min pre-gas + 2 min
PCU LED N ₂ (with nitrogen)	3 min at 50 % light output

Disposal

Disposal of contents/container in accordance with local/ regional/national/international regulations and according to the safety data sheet.

Hazards identification

Hazard and safety notes can be found on the product label and the corresponding safety data sheet.

Cleaning

The lacquered BTE earmolds and ITE shells should be cleaned before first use in patients. We recommend the OtoVita® care products series for this. Please observe the separate OtoVita® instructions. The use of aggressive cleaning agents can damage the material and the surface.

Disinfection

We recommend the OtoVita® care products series for disinfection of the lacquered BTE earmold or ITE shell. Please observe the separate OtoVita® instructions.

Sterilization

Sterilization of lacquered BTE earmolds or ITE shells is not required. NanoScreen™ Lack cannot be sterilized.

Serious incidents

All serious incidents occurring in relation to the product shall be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is settled.

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