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### Achtung!

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## EN Directions for use

### Safety advice!

The impression taking of the external auditory canal only has to be executed by qualified personnel. The working instructions and precautionary measures indicated in this working instruction have to be strictly observed. Non-observation might lead to irreparable damages of the ear or of the eardrum. Such still damages or not desired side effects appear despite the warnings, please send us a detailed description of your procedure and the damage or side effect at the patient.

### Intended use

Addition vulcanizing ear impression silicone.

### Indications for use

Ear impression material intended to reproduce a patients external auditory canal.

### Tarjet Patient Group

Patients who need ear care.

### Product description

Addition-vulcanizing ear impression material of soft, stable/thixotropic consistency for the slightly pressure-forming impression technique (Otoform® Ak, Otoform® Xpand) as well as for the impression technique without pressure (Otoform® A flex, Otoform® A soft, Otoform® A softX). While setting Otoform® there is a slight heat built-up. During ear impression with Otoform® Xpand the material property of mild expansion is illustrated by quietly crackling in the ear. Otoform® Xpand is not suited for deep impressions. Due to the platinum hardener system the material is nearly free of shrinkage and physiologically friendly.

### Processing

Before taking an impression the auditory canal and the eardrum have to be inspected. If any unusual features or abnormalities are found (e.g. irritation, perforation of the eardrum) an impression must not be taken. In cases of thick cerumen deposits and hairs in the auditory canal, the ear may require cleaning and the hairs trimmed. This may require referral to a suitable medical practitioner. In order to protect the eardrum, before the impression is taken an impression pad must be placed at the end of the external auditory canal in front of the eardrum. The material is automatically dosed out of the double cartridge in a mixing ratio of 1:1. Put the double cartridge into the injector. Remove the cap of the cartridge and dispense a bit of the material until both components flow equally. Now, attach the mixing canula onto the double cartridge. Make sure that the material is conveyed uniformly and mixed out of the mixing canula. For this, again a little material has to be conveyed before taking the actual ear impression. Slowly apply the impression material into the auditory canal until the complete concha is filled. Please observe that the mixing tip always touches the applied material. After use of the mixing canula leave the mixing canula on the double cartridge and do not exchange it before the next application. It serves as sealable lid. After the material has set the impression can carefully be removed from the ear. After the impression taking the auditory canal has to be inspected again. The Otoform® impressions can be disinfected with the OtoVita® Professional Disinfection Spray according to the instructions in the enclosed leaflet.

### Attention

The polymerization, i.e. the setting of addition-vulcanizing material can be interrupted by cerumen, latex gloves, cremes, resins etc. Therefore, we recommend to wear standard gloves made of polyethylene. Mixed material is chemically persistent. Spots onto clothing cannot be removed.

### Contraindication

In case of allergies to any of the ingredients do not use Otoform®. Not desired side effects of this medical device are not to be expected in case of proper processing and application. However, immune reactions (e.g. allergies) or irritations in the auditory canal principally cannot be excluded. Should you get to know non-desired side effects – also in doubtful cases – we appreciate in any case to be informed with a detailed description of the attendant circumstances and symptoms. We will follow each hint.

### Serious incidents

All serious incidents occurring in relation to the product shall be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is settled.

### Single-use advice

After impression taking, the cured silicone cannot be used for a further impression due to its physical features.

### Disposal

Disposal of contents / container in accordance with local / regional / national / international regulations and safety data sheet.

### Storage

Keep at room temperature 18°C–28°C (64°F–82°F).

### Lot number / Expiry date

The lot number and the expiry date are indicated onto the external packaging as well as on each double cartridge. In case of claims please always indicate the lot number of the product. Do not use the product the expiry date.

### For technical data please refer to the last page.

Dreve Otoplastik GmbH is not liable for any damages caused by improper application of the impression material.

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