Dubliermasse Duplicating material

6 kg

REF 40900

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Special indications for use: Too high melting temperatures can cause the duplicating material to darken, thus reducing the permeability to light and therefore diminishing the quality of polymerizing light-curing materials.

1. For use with duplicating units:

- Cut the duplicating material into small pieces and fill it into the duplicating unit EDG 1 (REF 61760 ff).
- Switch on the duplicating unit (refer to the operating instructions of the duplicating unit).

ATTENTION: Check the selected melting and processing temperature!

After filling the material into the investment form, several cooling processes can be used (see below).

2. For heating without duplicating unit EDG 1 (e. g. microwave oven, wax immersion and duplicating device, etc.):

During the heating process, **the temperature has to be checked continuously** with a thermometer.


- Cut the duplicating material into small pieces and heat to melting temperature.
- After reaching the melting temperature, the duplicating material has to cool to processing temperature. After filling the material into the investment form, several cooling processes can be used (see below).

If used according to the instructions, the duplicating material can be used 2x, depending on the heating method. Negative molds which are no longer needed, can be stored in the tightly closed original container until they are reused. They may not be put into the liquid duplicating material.

Cooling the negative molds:

1. In the refrigerator approx. 25 min.
2. In the freezer approx. 15 min.

Stability / Batch: See product.

Aybühlweg 59
87439 Kempten
Germany

Tel. +49 (0) 831 5 8113-0
Fax +49 (0) 831 5 8113-13

www.egger-labor.com