Zubehör / Accessories / Accessoires / Accesarios / Accessoires / Acessórios / Принадлежности / Akcesoria / 配套工具

1 x		REF 532, 533, 534
1 x		REF 4610
100 x		REF 5671, 5672, 5673, 5674, 5675
100 x		REF 565
100 x		REF 566

Technische Daten / Technical data / Données techniques / Datos técnicos / Technische gegevens / Dati tecnici / Dados técnicos / Технические данные / Dane techniczne / 技术数据

	Otoform® Ak	Otoform® Ak X
Verarbeitungszeit / Processing time / Temps de traitement / Tiempo de elaboración / Verwerkstijd / Tempo di lavorazione / Tempo de processamento / Время переработки / Czas przygotowania / 操作时间	100 sec ± 15 sec	100 sec ± 15 sec
Abbindezeit / Setting time / Temps de prise / Tiempo de endurecimiento / Verhardingstijd / Tempo di indurimento / Tempo di cura / Время склеивания / Czas wiązania / 硬化时间	3 min ± 15 sec	3 min ± 15 sec
Endhärte / Final hardness / Dureza final / Eindhardenheid / Dureza finale / Dureza final / Окончательная твердость / твердость конкава / 最终硬度	35 ± 2 Shore A	30 ± 3 Shore A
Verformung unter Druck / Strain in compression / Déformation sous pression / Deformación bajo presión / Verformung onder druk / Deformatie sotto pressione / Deformação sob pressão / Изменение формы под давлением / Однородность под давлением / 压力下的变形	4.3 % ± 0.5 %	4.4 % ± 0.5 %
Rückstellung nach Verformung / Elastic recovery / Recuperación tras deformación / Reposition tras deformación / Regeneracja po odkształceniu / 变形后复位	> 99.6 %	> 99.6 %
Lineare Maßänderung / Linear dimensional change / Variation dimensionnelle linéaire / Cambio dimensional lineal / Lineaire maatverandering / Modifica dimensionale linea / Mudança dimensional linear / Линейное изменение размера / Liniowa zmiana wymiaru / 测量的线性变化	< 0.2 %	< 0.2 %
Konsistenz / Consistency / Consistenza / Consistencia / Consistenza / Consistência / Консистенция / Konsystenca / 稳定性	knetbar – Typ 0 putty – type 0 malleable – type 0 viscosidad elevada – tipo 0 viscoelasticidad elevada – tipo 0 impastabile – tipo 0 calafatado – tipo 0 plástica – tipo 0 daający się ugniatać typ 0 柔软可揉捏 类型 0	knetbar – Typ 0 putty – type 0 malleable – type 0 viscosidad elevada – tipo 0 viscoelasticidad elevada – tipo 0 impastabile – tipo 0 calafatado – tipo 0 plástica – tipo 0 daający się ugniatać typ 0 柔软可揉捏 类型 0

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EN

Caution!

The impression of the outer ear canal may only be taken with impression material by persons with appropriate qualifications. The processing information and precautions which are set out in these directions for use must be observed under all circumstances. Failure to observe these may result in irreparable damage to the ear or the eardrum. If damage or undesirable side effects occur despite the warnings, we request that you report this with a detailed description of your procedure and the damage to or side effects for the patient.

Intended use

Additional curing ear impression silicone

Indication for use

Ear impression material for reproduction of the outer ear canal of a patient with hearing loss, Tinnitus aurium or impaired hearing for the provision and adjustment of a hearing device.

Contraindication

The product may not be used if the patient is allergic to one of the ingredients. There are no side-effects to be expected from this medical device if processed and used properly. However, immune reactions (e.g. allergies) or local disorders of sensation (e.g. irritation in the auditory canal) cannot be ruled out in principle. If you

become aware of adverse side-effects, also in cases of doubt, please be sure to advise us in any case and describe the circumstances and symptoms as accurately as possible. We will investigate any notice.

Patient groups

Otoform® must not be used for patients suffering from a perforated eardrum or injuries to the ear canal.

Target Patient Group

Patients who need ear care.

Intended users

ENT specialist, hearing aid audiologist

Product description

The ear impression material consists of addition curing silicone with a soft, thixotropic consistency for impression technology which creates slight pressure. The material is practically shrink-proof and is physiologically well tolerated thanks to the platinum hardener system. Slight generation of heat can be felt while Otoform® is curing.

Processing

The ear canal and the eardrum must be examined before the impression is taken. No impression may be taken in the event of abnormal findings such as inflammation or perforation of the eardrum. Cleaning of the ear canal and removal of hairs by an ENT doctor is required in the event of significant earwax deposits and hairs in the ear canal. An impression pad must be placed at the end of the outer ear canal before the eardrum before the impression is taken in order to protect the eardrum. The two components of the impression material are removed from the containers in equal parts with measuring spoons and pressed onto a mixing tray. The two components are now rapidly kneaded together with a spatula until the mixture is a uniform colour. For hygiene reasons, kneading by hand is only recommended with gloves (no latex gloves!). The kneaded mass is now filled into an impression syringe. The impression material is slowly applied into the ear canal and the concha. Please take care to ensure that the tip of the impression syringe is always in direct contact with the impression mass during application. Once the material has hardened, the impression can be carefully removed from the ear. The ear canal must be examined again after the impression has been taken. Carefully close the tips after use.

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Disposal

Disposal of container in accordance with local / regional / national / international regulations and safety data sheet.

Hazard identification

Hazard and safety notes can be found on the product label and the corresponding safety data sheet.

Serious incidents

All serious incidents occurring in relation to the product shall be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is settled.

Batch number / use-by date

The batch number and the use-by date are located on both the outer packaging as well as on each tub. Both components A and B of a packaging must only be used in this combination. When submitting a complaint about a product, please always quote the batch number. Do not use the product after expiration of the use-by date.

See the last page for technical data and required specifications.

Dreve Otoplastik GmbH is not liable for any damages caused by improper application of the product.

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