egger A/soft egger A/ITE egger A/II

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General information

egger A/soft

Silicone impression material with an extremely easy flow and a high degree of form stability, non-sticky and smooth. Color: white/blue.
Final degree of hardness: 35 shore A

egger A/ITE

Silicone impression material with a very soft consistency, an extremely easy flow and a high degree of form stability. Impression taking with minimum pressure possible, especially suited for CIC impressions and those of small ear canals. Color: white/pink.
Final degree of hardness: 34 shore A

egger A/II

Silicone impression material of slightly pressure-forming consistency. Due to its high viscosity, the material is characterized by impression exactitude, tensile strength and form stability. Color: red.
Final degree of hardness: 40 shore A

Due to a high-quality platinum hardening system, the materials are virtually shrink-proof and are characterized by their highly physiological compatibility.

| | | |
|-----------------|--------------|----------------|
| | | low _____ high |
| Flow- | egger A/soft | _____ |
| ability: | egger A/ITE | _____ |
| | egger A/II | _____ |

Packaging units for egger cartridges

Basically, there are 8 double cartridges per packaging unit. egger A/soft and egger A/ITE are available in both cartridge systems, whereas egger A/II is only available in the S 50 cartridge system, due to its high viscosity.

| | | |
|--------------|--------------------|-----------|
| egger A/soft | Standard cartridge | REF 20600 |
| | S 50 cartridge | REF 20617 |
| egger A/ITE | Standard cartridge | REF 20604 |
| | S 50 cartridge | REF 20608 |
| egger A/II | S 50 cartridge | REF 23707 |

Required accessories for standard cartridge systems

| | | | |
|-------------------------|-----------|--------------------------------|-----------|
| Inject pistol | REF 25100 | | |
| Mixing canula Ø 6.75 mm | | Mixing canula »Slim« Ø 4.75 mm | |
| 50 pieces | REF 61403 | 100 pieces | REF 61420 |
| 100 pieces | REF 61402 | | |

Required accessories for S 50 cartridge systems

| | | | |
|---|------------------|-----------------|--|
| Inject pistol | REF 25150 | | |
| Mixing canula S 50 for egger A/soft and egger A/ITE | | | |
| | yellow, Ø 5.7 mm | pink, Ø 6.75 mm | |
| 40 pieces | REF 25160 | REF 25151 | |
| 100 pieces | REF 25161 | REF 25152 | |
| Mixing canula S 50 for egger A/II | | | |
| green, Ø 8 mm | 40 pieces | REF 25154 | |
| | 100 pieces | REF 25155 | |

Furthermore, we can offer you a wide range of diagnostic instruments such as otoscopes and ear lights as well as impression pads in different sizes and designs.

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






Product information:

Storage: The material is to be stored in a dry and cool place. Keep the cartridges always closed.

Stability / Batch:

See packaging and/or on the cartridge.

Mixing: The two components are mixed automatically by the mixing canula 1:1. Then, the material immediately is ready for processing.

Curing times (apply for all materials): 4–5 min.

The processing times refer to a room temperature of $23\text{ °C} (73.4\text{ °F}) \pm 2\text{ °C} (35.6\text{ °F})$ and a relative air humidity of 50%, the hardening time refers to a body temperature of $37\text{ °C} (98.6\text{ °F}) \pm 2\text{ °C} (35.6\text{ °F})$. Lower temperatures prolong, higher temperatures decrease the processing and hardening times.

Processing

1. Examining the ear: By making the otoscopy, anomalies in the ear and/or in the auditory canal area can be seen. If distinctive features are found, do not take any impressions and seek medical advice.
2. Preparing the ear: Depilate the auditory canal, clean the impression area if necessary. Position the impression pad properly to protect the eardrum.
3. Insert double cartridge into inject pistol, take off sealing cap of the cartridge, then take out a pea-sized quantity and attach the mixing canula. Creaming the auditory canal (e.g. baby care oil) facilitates demolding.
4. Impression taking: Please take care that the injector is safely supported.
5. Demolding: After the specified time, carefully take out the impression. Please mind the pressure compensation. The length of the impression should include the second bend of the auditory canal part.
6. Check: Check the ear and auditory canal after the impression taking with the otoscope.

Notes:

- To guarantee a constant mixing ratio, use the inject pistol to form a pea-sized amount. Then, attach the mixing canula.
- In case the inject process is interrupted for a longer period of time, the mixing canula has to be replaced.
- Already used cartridges can be stored and used until they are empty.
- Do not use latex gloves. Latex contaminated surfaces, cerumen, creams and resins can influence the hardening. The two components may not be mingled before mixing them.

Caution: Cross-linked materials are chemically-resistant – avoid stains on clothes!

Important: The impression taking is to be carried out by qualified personnel only! Please adhere to the processing instructions and safety regulations that are described in the product information.