**Scopo d’uso**

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**Intended use**

Light-curing material for the production of earmolds.

**Extended purpose**

Light-curing material as an accessory for the production of earmolds.

**Indication**

Hearing loss or impaired hearing

**Indication for use**

Polylux 1000N Gel for repairs or modelling of vents for BTE earmolds and ITE shells.

**Contraindication**

If an acrylic allergy exists, an allergic reaction may occur.

**Patient target group**

People with impaired hearing ability to wear a hearing aid.

**Intended users**

Qualified personnel from earmold laboratories and hearing aid acousticians.

**Requirements / equipment**

Polylux 1000N

**Basic Material Properties**

See technical datasheet (TDS). Available on request.

**Specific Manufacturing Considerations****Caution**

During the application we recommend wearing personal protective equipment i.e. suitable gloves, safety goggles etc. Use a suitable extraction system or ensure sufficient ventilation when working with the lacquer.

**Further processing****Preparation**

The surfaces of the BTE earmold or the ITE shell should be free from dust and grease residues. Roughening of the plastic surface is recommended for improved wetting behaviour. Installing the dosing syringe

**Lot number / Expiry date**

The batch number and the use-by date are located on both the outer packaging and the product. When submitting a complaint

**Further processing**

Fotoplast® Gel is suitable for repairs, for modelling of a venting or similar work on hard earmolds.

Repair: Apply Fotoplast® Gel and allow to cure (➔ see post-curing).

Venting: Place venting wire and cover with Fotoplast® Gel, then allow to cure (➔ see post-curing) and remove venting wire.

**Post-curing**

Polylux 1000N (with nitrogen) clear, red and blue-transparent	2 min
PCU LED N <sub>2</sub> (with nitrogen) opaque H (yellowish)	6 min at 80% light output

**Disposal**

Disposal of contents/container in accordance with local/ regional/national/international regulations and according to the safety data sheet.

**Hazards identification**

Hazard and safety notes can be found on the product label and the corresponding safety data sheet.

**Cleaning**

The BTE earmolds and ITE shells should be cleaned before first use in patients. We recommend the OtoVita® care products series for this. Please observe the separate OtoVita® instructions. The use of aggressive cleaning agents can damage the material and the surface.

**Disinfection**

We recommend the OtoVita® care products series for disinfection. Please observe the separate OtoVita® instructions.

**Sterilization**

Fotoplast® Gel cannot be sterilized.

**Serious incidents**

All serious incidents occurring in relation to the product shall be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is settled.

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