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Temperaturregler TR4

REF 66250

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
REF 66250 (230 V)

REF 66251 (115 V)

Please read the operating instructions carefully before working with the unit! Please keep the operating instructions in a safe place for future reference!

	<p>Important Information</p>
	<p>Caution</p>

<u>Index</u>	Page
1 General information	17
2 Safety regulations	17
2.1 Working according to the safety regulations	17
2.2 Usage according to the regulations	18
3 Installation	18
3.1 Scope of delivery	18
3.2 Technical data	18
3.3 Spare parts and accessories	18
3.4 Assembling and connecting the unit	18
4 Operation	19
4.1 Operational and control elements	19
4.2 Selecting the heating output	19
4.3 Data storage	20
4.4 Parametric adjustment of the desired temperature values	20
4.5 Adjustment of the set value	20
4.6 Description of the adjustment procedure	21
4.7 How to exit the parametric level	22
4.8 List and explanation of the parameters	22
4.9 Calibration of the sensor	26
5 Putting into operation	26
5.1 Waxing of ear impressions	26
6 Maintenance	27
6.1 Refilling and exchanging the immersion wax	27
7 Warranty/Service/After-sales service	28

1 General information

The temperature control unit TR4 is designed for monitoring the temperature of the wax. It is inserted between the wax melting appliance and the socket in order to guarantee a constantly controlled temperature.

The heating power of the unit can be individually adjusted by means of a potentiometer. The temperature sensor which is connected to the temperature control unit is located inside the wax. It controls the heating power of the wax melting unit.

The temperature control unit TR4 allows entering two different desired temperature set values. After turning on the unit, the control system heats up to the indicated set value. Once this value is achieved, the control unit switches over from a 100% heating power to the heating output adjusted by means of a potentiometer.

The temperature control unit was built according to state-of-the-art technology and is electrically protected.

The appliance is easy to use and maintain.
The elegant solid metal casing is easy to clean.

2 Safety regulations

	<p>The temperature control unit TR4 is used in conjunction with a wax heating unit which was also built according to state-of-the-art-technology. A usage other than described in the operating instructions may lead to injuries of the operator or of third persons. Furthermore, the device or other material assets may be damaged. Therefore, the device must be used in technically perfect condition only and the safety regulations must be strictly complied with.</p>
---	---

2.1 Working according to the safety regulations

	<ul style="list-style-type: none"> • The device must be used by qualified and trained persons only (especially children must be kept away from the unit!) • Before putting the device into operation, the user has to read the operating instructions and ask the manufacturer in case of questions. • Connect the device to permissible sockets only, see technical data, section 3.2. • Before turning on the unit, the sensor absolutely has to be immersed in the wax. • Use the unit on dry working places only. • Use the unit in clean condition only. The current instructions for the material to be worked on must be strictly observed. (Otherwise quality divergence may occur). • The device must not be opened by unauthorized persons. Before opening the device, pull the mains plug. • Protect the unit from dust. • Should the mains cable be damaged, it must be replaced at once. Defective mains cables must not be used. • Do not pull the mains plug out of the socket at the mains cable. • The device must be repaired by qualified persons only. • The unit must be switched off during service and maintenance operations (cleaning, exchange of the sensor or the dipping wax).
---	---

2.2 Usage according to the regulations

Only waxes may be heated up in the wax heating unit/control unit TR4.
This combination of appliances is designed for waxing ear impressions only.

Do not use any dissolvent when cleaning the device. We generally recommend using surface cleaner (REF 31601 ff). Do not immerse the unit in liquids.

- This combination of devices must only be used for waxing ear impressions.

In case of self-made repairs or use of other spare parts all guarantee and warranty claims will become extinct!

3 Installation

3.1 Scope of delivery

Please check whether the delivery is complete and whether all parts are free from transport damage. Immediately inform the distributor about missing parts or transport damage.

1 control unit with temperature sensor
1 operating instructions

3.2 Technical data

Power supply: 230 V +/- 10% (or 115 V)
Power frequency: AC 50 Hz
Connection: German shock-proof plug
Power input: 2.3 kW

3.3 Spare parts and accessories

Temperature sensor
Operating instructions

3.4 Assembling and connecting the unit

The temperature control unit corresponds to the state-of-the-art standards and to the acknowledged safety regulations, nevertheless dangers can arise. Please observe the following before putting the unit into operation:

- The unit must be operated by qualified persons only.
- Observe the recognized industrial safety regulations.
- The unit must be used in technically perfect condition only and in compliance with the operating instructions.
- Position the device horizontally and place it in its suitable working position.
- Protect the unit against impact.
- Check whether the unit is dry.

- Check whether the device is in technically perfect condition, especially the mains cable!
- Connect the temperature sensor to the screw socket at the front panel of the device.
- Put the temperature sensor into the wax heating unit.
- Plug the mains cable of the wax heating unit into the socket that is located on the rear side.
- Please bear in mind that only consumers with a maximum output of 3300 W can be connected.
- Plug in the mains cable.

4 Operation

The temperature control unit TR4 is used in connection with a wax jug to heat immersion wax.

4.1 Operational and control elements

Main switch with green LED

Position: Front panel of the device.

Function: Power supply, should be turned off after using the device.

Potentiometer

Position: Front panel of the device.

Function: It controls the percentage of heating capacity of the wax heating unit (1 – 99%).

Screw plug box for temperature sensor

Position: Front panel of the device.

Function: Connection of the temperature sensor.

Display with "SET" and "↑ ↓" keys

Location: Frontal panel of the device.

Function: If you simultaneously press the two "↑ ↓" keys for approx. 3 seconds, the control unit automatically switches over to the parametric level (P1 - P32). If no button is pressed, the display shows which parameter of the list is selected, starting with P1. To switch over to different parameters, press the "↑" or "↓" key. Pressing the "SET" key, the desired parameter appears in the display. Pressing the "SET" key and the "↑" or "↓" keys simultaneously, you can increase or reduce the parameter. As soon as you let go the "↑" or "↓" key, the value is stored. If you want to switch from the parametric level over to the set value level, press the "↑" and the "↓" key at the same time for 3 seconds, or just wait 60 seconds.

4.2 Selecting the heating output

The control unit is provided with an electronic potentiometer which allows adjusting the heat very fast. (The power is controlled by a full wave control unit). The load resistance (heating) is switched ON and OFF at intervals. Modifying the "ON" and "OFF" cycles via the turning knob you can adjust the operating time (1 – 99 %). If the operating time is perfectly set, a narrow adjustment range is the result. The resistance load is switched free of interferences (zero crossing). The yellow LED (cycle) shows the function of the potentiometer and flashes.

Particular feature: Parameter P1 (see adjusting the parameters) corresponds to a second set value (set value 2). The unit heats up until this set temperature is achieved using 100% of its capacity and then switches over to the selected heating power / potentiometer cycle.

<u>Example:</u>	Set value 1	95 °C (203 °F)
	Potentiometer	25%
	Set value 2	80 °C (176 °F)

The desired final temperature of the material is 95 °C (203 °F).

The device heats up the material with 100% of its heating capacity until it achieves the temperature of 80 °C (176 °F). If the temperature of 80 °C (176 °F) is reached, the electronic system reduces the heating power to 25% - the result is a small adjustment range.

4.3 Data storage

All parameters that are entered via the keyboard, are saved to an EPROM and are safeguarded even in case of a sudden power outage.

4.4 Parametric adjustment of the desired temperature

The control unit can be adjusted (parameterized) on two different operating levels. The first level is called the set value level and allows the access to all adjusted parameters; on the second level, which is called the parameter level, all parameters can be adjusted. We deliberately complicated the access to the parameter level in order to avoid an erroneous setting of the parameters.

4.5 Setting the set value 1

Display of the currently measured temperature

Adjustment of the desired wax temperature

It can be adjusted within the nominal value range set (see P4 and P5) press the "SET" key and the desired value can be adjusted by pressing the "↑" or "↓" keys. Factory-made, the wax temperature is set to 95 °C (= 203 °F).

Before putting the temperature control unit TR 4 into operation for the first time, the following procedures are necessary:

- ① Calibration of the sensor
- ② Parameterisation P1-P32
- ③ Adjustment of the set value 1

Please note that

- the electric wax jug (REF 65005) is set to its maximum heating power (= setting 6).
- the electric wax jug is continuously filled with at least $\frac{3}{4}$ of green immersion wax (REF 41201).
- the user instructions only apply with the use of the egger immersion wax (REF 41201).
- the use of other waxes may cause changes in the parameters
- the already factory-made settings, under normal circumstances, guarantee an immediate use. In case you would like to change the parameters, please adhere to the following notes.

4.6 Description of the adjustment procedure

First, you have to press the keys "↑" and "↓" at the same time for approx. 3 seconds. Then, the control unit switches over to the parameter level with the parameter list P1 - 32. If no key is pressed, the display shows the selected parameter of the list, starting with P1. If you want to change over to different parameters, just press the "↑" or "↓" key. By pressing the "SET" key, the adjusted parameter appears on the display.

If you simultaneously press the "SET" key and the "↑" or "↓" key, you can increase or reduce the parameter. Once you let go the "↑" or "↓" key the parameter is automatically stored.

4.7 How to exit from the parameter level

If you want to go back to the set value level, simultaneously press the “↑” and the “↓” key for 3 seconds or just wait 60 seconds, then the control unit automatically switches over to the set value level.

4.8 List and description of the parameters

With the “P” settings, the basic parameters can be changed. To get to level “P” you have to press simultaneously the “↑” and the “↓” key for about 2 seconds.

- P1 Setting of the full load temperature
- P2 Adjusting the limit of tolerance downwards
- P3 Adjusting the limit of tolerance upwards
- P4 Adjusting the nominal value downwards
- P5 Adjusting the nominal value upwards
- P6 Calibration of the temperature sensor
- P19 Key interlock
- P30 Adjusting the optical limit value alarm if underran
- P31 Adjusting the optical limit value alarm if exceeding
- P32 Adjusting the limit of tolerance for the limit value alarm

“P1” = Setting of the full load temperature

Up to this temperature (for example 85 °C/185 °F) the temperature control unit heats up (to 100% heating power) under full load. As soon as the set temperature is reached, the device reduces the heat pulsing. Factory-made, the temperature is set to 80 °C(176 °F).


We recommend adjusting the temperature (“heating power in %”) to 25% by using the red switch!


“P2” = Adjusting the limit of tolerance downwards

With this option, it is possible to define downwards the maximum deviation value (less than the nominal value). In the example shown, the maximum deviation with 0.3 °C (32.54 °F), which means that at a nominal value of 95 °C (203 °F) a maximum lowest value of 94.7 °C (202.46 °F) is accepted by the user.


“P3” = Adjusting the limit of tolerance upwards

The same function as described for “P2”. Here, however, the maximum deviation value is defined upwards (greater than the nominal value). In the example shown, the maximum deviation of 0.3 °C (32.54 °F) is adjusted, which means that at a nominal value of 95 °C (203 °F) a maximum peak value of 95.3 °C (203.54 °F) is accepted.


“P4” = Adjusting the nominal value downwards

With “P4” the maximum value of the wax temperature is defined downwards. Hence, the wax temperature, if required, can be changed downwards fast at the push of a button as described under 4.9. Factory-made, the temperature is adjusted to 90 °C (194 °F).


“P5” = Adjusting the nominal value upwards

As done with “P4” downwards, here, with “P5” the maximum value of the wax temperature is defined upwards. Hence, as described under 4.9, the wax temperature can be changed upwards fast at a push of a button. Factory-made, the temperature is set to 95 °C (203 °F).


“P6” = Calibration of the temperature sensor

With this option, differences in the temperature sensor can be compensated (caused by sea level air pressure). Factory-made, the value is reset to zero.


“P19” = Key interlock

With the setting to “1” the key interlock is activated. This avoids an unintentional change of the parameters. By setting it to “0” the key interlock is deactivated.


“P30” = Adjusting the optical limit value alarm

With this setting you can activate an optical alarm (blinking) as soon as the desired wax temperature is under-run by the value set. This example shows how the optical alarm is triggered as soon as the nominal value is under-run 2 °C/35.6 °F, for example at 93 °C (199.4 °F) if the nominal value is 95 °C (203 °F).

“P31” = Adjusting the optical limit value alarm

As described at P30, but the activation of the optical alarm is effected by exceeding the nominal value. This example shows how the optical alarm is triggered as soon as the nominal value is under-run by 5 °C/41 °F, for example at 100 °C (212 °F), if the nominal value is 95 °C (203 °F).

“P32” = Adjusting the limit of tolerance for the limit value alarm

This function is not applicable. Factory-made, the value is set to 0.3.

4.9 Calibration of the sensor

The temperature sensor is calibrated on the occasion of the initial operation or after having replaced a faulty sensor. Boil some water in a pot and dip the temperature sensor in it. The display now shows the actually measured temperature (set 1) of the water, which should be 100 °C (212 °F).

If the sensor indicates a different temperature, change to the programming level P6 (see the parameterization chapter) where you can enter the corresponding correction factor in order to compensate the difference. Then, go back to the set value level and remeasure the temperature of the water. If there are still deviations from the desired 100 °C (212 °F), you can switch over again to the parameterization level P6 in order to correct the value as often as you want.

	Functional description	Setting range	Wax heating unit
P1	set value 2		80 °C (176 °F)
P2	hysteresis (dual-mode)	0,1... 99K	0.3
P3	hysteresis (dual-mode)	0,1...99K	0.3
P4	minimum set value	-99...999 °C (-146.2...1830.2 °F)	90 °C (194 °F)
P5	maximum set value	-99...999 °C (-146.2...1830.2 °F)	110 °C (230 °F)
P6	actual value correction	-10...+10K	depending on the sensor
P19	keyboard locking	0: not locked 1: locked	1
P30	minimum limit alarm		- 2
P31	maximum limit alarm		+ 5
P32	hysteresis limit alarm		0.3

5 Putting into operation

	<p>The appliance is exclusively designed for the waxing of ear impressions. Never operate the device without immersion wax!</p> <p>The temperature control of the wax heating unit always has to be set to maximum output, the power supply is adjusted via the temperature control unit TR4. The temperature sensor must always remain in the wax heating device!!</p> <p>Please bear in mind that failure may arise despite the state-of-the-art technology used. Please ensure that an overheating protection such as a temperature controller, temperature limit or temperature fuse will come into force in case of a possible failure of the temperature control unit; or operate the device only under constant supervision. The manufacturer does not accept any liability for resulting damages.</p>
---	---

5.1 Waxing of ear impressions

The following conditions must be fulfilled:

The mains plug must be plugged in.

Switch on the device (main switch).

Please bear in mind that the unit takes a certain time to heat the wax (the yellow lamp is illuminated).

Before starting the waxing of the impressions, please check the wax temperature.

Put the impression on the dipping sieve and immerse it steadily in the wax for a short time and then take it out again.

After a short phase of cooling down, the waxed impressions can be taken off the dipping sieve.

	<p>Please note that</p> <ul style="list-style-type: none"> - the electric wax jug (REF 65005) is set to its maximum heating power (= setting 6). - the electric wax jug is continuously filled with at least $\frac{3}{4}$ of green immersion wax (REF 41201). - the user instructions only apply with the use of egger immersion wax (REF 41201). - the use of other waxes may cause changes in the parameters - the factory-made settings, under normal circumstances, guarantee immediate use. In case you would like to change the parameters, please adhere to the notes for parametric adjustments.
---	---

6 Maintenance

	<p>Before every maintenance work, the device must be separated from the mains supply. Hot oil – danger of burning! Move or transport the unit only if it has cooled down.</p>
--	---

6.1 Refilling or exchanging the immersion wax

If you regularly wax ear impressions, the wax level in the wax container gradually drops. Be sure to have a sufficient amount of wax on stock.

Fill the pieces of the wax into the wax heating container.

It takes some time for the wax to heat to the correct temperature.

The unit is now ready to be used again and can be connected to the power supply.

7 Warranty/Service/After-sales service

Warranty

The warranty period for our devices starts with the date of delivery and is 12 months. Excluded from the warranty are modifications of the unit and all other parts of the unit which are subject to wear and tear.

	In case of improper handling of the temperature control unit TR 4 all guarantee and warranty claims will become extinct!
--	--

We will be pleased to be of assistance to you in case of damage or other malfunctions:

egger Otoplastik + Labortechnik GmbH
Aybühlweg 59
87439 Kempten
Phone: +49 831 58113-60
Fax: +49 831 58113-14
E-mail: labortechnik@egger-labor.de
Internet: www.egger-labor.com

Please indicate the following data:

Device type:
Temperature control unit TR 4
(REF 66250/230 V or 66251/115 V)

..... (see label)

Serial-no.:

..... (see label)