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Achtung!

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Directions for use

Safety advice!

The impression taking of the external auditory canal only has to be executed by qualified personnel. The working instructions and precautionary measures indicated in this working instruction have to be strictly observed. Non-observance might lead to irreparable damages of the ear or of the ear drum. Should still damages or not desired side effects appear, please send us a detailed description of your procedure and the damage or side effect at the patient.

Product description

The ear impression material consists of a condensation-vulcanising silicone with a soft, non-sticky consistency and is especially suitable for functional impression-taking. During the impression taking the material exerts a slight pressure onto the ear tissue being advantageous for the fitting and the sealing of the earmold.

Total processing time:	≥ 65 sec
Setting time:	4 min ± 30 sek
Final Hardness:	30 ± 2 Shore A
Deformation under pressure:	5.8 ± 0.5 %
Recovery from deformation:	> 99.8 %
Linear dimensional change:	< 0.3 %
Consistency:	putty - Type 0

! The mixing and processing times refer to a room temperature of 23 ± 1 °C and a relative air humidity of 50 %, the setting time to a body temperature of 37°C ± 1°C.

Processing

Before taking an impression the auditory canal and the ear drum have to be inspected. If suspicious details are found e. g. irritations, perforations of the ear drum an impression must not be taken. In case of thick cerumen deposits and hairs in the auditory canal it has to be cleaned and the hairs have to be eliminated. In order to protect the ear drum an impression pad has to be placed at the end of the external auditory canal in front of the ear drum before an impression is taken. Please take basis material out of the tub with the help of the dosing spoon. The dosing spoon must be filled with basis material up to the top. Put the material onto a mixing block and form it to a plate. Press the dosing spoon onto the formed plate in order to get the dimension for the hardener paste. Press a line of approx. 4.5 cm of the hardener paste out of the tube (diameter of the tube opening) onto the formed plate of the basis material. The long side of the spoon impression corresponds to 4.5 cm. Quickly mix both components with a spatula until a homogeneous colour is achieved. Fill the mixture into an impression syringe and apply the mass into the auditory canal and the concha. Please observe that the mixing tip always touches the applied material. When the material has set, the impression can carefully be taken out of the ear. After the impression taking the auditory canal has to be inspected again.

Attention

By overdosing of the hardener paste the setting time is accelerated and the shrinkage of the material is increased.

Hardener paste

Contains: Dioctyltinacetylacetonate, tetrakis (2-butoxyethoxy) silane, dioctyltin oxide, trimethoxyvinylsilane
Before use observe the Material Safety Data Sheet of Hardener paste.

Hazard statements

H315 Causes skin irritation. H317 May cause an allergic skin reaction. H371 May cause damage to organs. H373 May cause damage to organs through prolonged or repeated exposure.

Precautionary statements

P260 Do not breathe dust/fume/gas/mist/vapours/spray. P280 Wear protective gloves. P362+P364 Take off contaminated clothing and wash it before reuse. P333+P313 If skin irritation or rash occurs: Get medical advice/attention. P405 Store locked up. P501 Dispose of contents/container in accordance with local/regional/national/international regulations.

Warning

Side effects

The direct skin contact with the hardener paste has to be avoided as the paste is irritating.

Not desired side effects of this medical device are not to be expected in case of proper processing and application. However, immune reactions (e. g. allergies) or e. g. irritations in the auditory canal principally cannot be excluded. Should you get to know non-desired side effects – also in doubtful cases – we appreciate in any case to be informed with a detailed description of the attendant circumstances and symptoms. We will follow each hint.

Form of delivery

REF	500	124 ml tub with 5 g hardener
REF	502	584 ml tub with 20 g hardener incl. dosing spoon
REF	504	13,1 l bucket without hardener
REF	458	1 dosing spoon, red

Storage

Keep well closed at room temperature.

Lot number/Best before date

The lot number and the best before date are indicated onto the external packaging as well as on each tub and bucket. In case of claims please always indicate the lot number of the product. Do not use the product after expiry of the best before date.

Dreve Otoplastik GmbH is not liable for any damages caused by improper application of the impression material.