

Hardener Paste



Technical Data Sheet

Intended use: Condensation curing ear impression silicone.

Manufacturer: Dreve Otoplastik GmbH
Max-Planck-Straße 31
59423 Unna / Germany

Composition:

Pos	chemical term
1	organic tin compounds
2	silicic acid ester
3	vinyltrimethoxysilane
4	food colour
5	silicon dioxide
6	vaseline

Form: paste

Colour: red
(corresponds to released sample in the house)

Density: 1.02 g/cm³ at 20 °C DIN 51757

Solubility in water: insoluble

pH-value: neutral

Dreve Otoplastik GmbH

Max-Planck-Str. 31

59423 Unna / Germany

Tel.: 02303 / 8807-0

Fax: 02303 / 82909

otoplastik@dreve.de

www.dreve.com