

Thermosoft Lack



Technical Data Sheet

Intended Use: Light-curing lacquer for earmolds

Manufacturer: Dreve Otoplastik GmbH
Max-Planck-Straße 31
59423 Unna / Germany

Composition:

Pos	chemical term
1	acrylates
2	methacrylates
3	urethan(meth)acrylates
4	photoinitiator mix
5	additives

Color: clear transparent

Density at 20 °C: approx. 1.03 g/cm³

Viscosity at 23 °C: approx. 0.03 Pas

Dreve Otoplastik GmbH

Max-Planck-Str. 31

59423 Unna / Germany

Tel.: 02303 / 8807-0

Fax: 02303 / 82909

otoplastik@dreve.de

www.dreve.com